

ChiroSolve

**Enantiomer preparation for discovery,
development and manufacturing**

Presenter: Dr. Niteen A. Vaidya

CTO, ChiroSolve Inc.

2065 Martin Ave. suite 107,

Santa Clara, CA 95050

Telephone: (408) 834-8597

Fax: (408) 351-7900

Email: neelamv@chirosolve.com

URL: <http://www.chirosolve.com>

Background on Chiral Chemistry

Chirality in pharmaceutical market


- Chiral molecules are closely related, mirror-image molecules with +ve or –ve optical rotation; very often with different bioactivity
- Worldwide single enantiomer drugs market is over \$200 billion; 40% of total drug market
- 9 out of 10 top selling drugs are chiral in nature
- Number of chiral molecules as new drug candidates entering pre-clinical stage per year: 400
- Number of chiral drug molecules being considered as potential candidates: 40,000

Why chirally pure drugs?

Industry need to	How chirally pure product solves this need
Ensure product* safety	Unwanted enantiomer can be toxic, may cause side-effects
Ensure product* efficacy	Enantiomers have different (sometimes opposite) bioactivity
Minimize dosing	Enantiomerically pure products cut dosage by up to 50%, increase patient tolerance
Save time and money	Pure enantiomer saves clinical testing time, reduces bring-to-market cost – no requirement to test racemate or the other enantiomer
Comply with regulations	FDA and EU require fully study and disclosure of bio-activities of each enantiomer in drug; proof of enantiomeric purity using prescribed methods

Pathways to pure enantiomers

Chiral technologies



65% of chirally pure products are manufactured using diastereomeric crystallization. Chirosolve products are based on Diastereomeric crystallization technique


Diastereomeric crystallization: example separation

(+/-)Ibuprofen


(-) phenethyl
amine

2-propanol


Diastereomeric salt
(-) crystalline and (+) filtrate




+




+





=




Treat with strong acid and separate out



Scale-up to identify
of re-crystallizations
needed to get target purity



(-)isomer (crystalline)



(+)isomer (filtrate)

Diastereomeric crystallization challenges

- Identifying the optimal separation using diastereomeric crystallization technique involves choosing the right reagent and solvent out of 100s of combinations.
- High cost of initial material during discovery stage limits the ability to explore chiral separation routes effectively
- Due to time and resource constraints, pre-selecting small set of reagents to screen the racemate leads to sub-optimized method development
- Manual approach to both initial separation and optimization is monotonous, time consuming, and error-prone

ChiroSolve Products

Solution: ChiroSolve Products

- **Screen**
 - Exhaustive screening of racemate against 384 combinations provides comprehensive separation conditions for IP coverage
 - Short-lists the ideal separation conditions to be explored during method development
- **EnantioPrep (Screen + Purification kits)**
 - Isolates pure enantiomer from racemate for discovery and preclinical
 - Produces small quantity (up to 20 gm) of enantiomer within matter of days
- **ScalePrep (screen + Purification + Optimization kits)**
 - Explores short-listed separation conditions to identify shorted route to produce pure enantiomer
 - Further optimizes selected condition to identify ideal ratio between racemate and reagent; ideal concentration of solvent(s)
 - Optimized method is scalable to kilo quantities

ChiroSolv[®] advantage

- Short Term
 - Exhaustive and parallel screening provides several suitable options
 - Consistent results; robust method
 - Considerable time saving, resource cost saving
 - High versatility increases project success
 - Well documented step-by-step procedure eliminates human errors
- Long Term
 - Minimize research and manufacturing costs by finding the optimum strategy to employ when dealing with chiral separation
 - Reduce time-to-market by streamlining your research processes
 - Low cost, robust method to follow
 - Expand research options
 - Unbiased results ensure full exploration of separation methods
 - Establish IP foundation

ChiroSolv[®] Screen kits

- **384 combinations of resolving agents/solvents to screen against in parallel**
 - 4 kits per racemate; each with 96 vials (4 acid kits, 4 base kits)
 - Consistent results; identifies comprehensive set of conditions
 - Exhaustive screening offers unbiased results
- **Very little racemate needed** per experiment (0.001 to 0.03 mmol per vial)
 - Researchers do not have lot of sample to work with during early research
- **High-throughput standard kits**, compatible with most auto-stations
 - Vials and rack are chemically inert and withstand extreme temperature
 - Whole experiment can be done inside the rack without taking out vials
 - No additional setup required; kits are ready to use out of the box

ChiroSolv[®] Screen kits (cont')

- **Designed to screen any racemate, giving results within 24 hours**
 - Racemate can be acid, base, alcohol, aldehyde, ketone or amino acid and can have single or multiple chiral centres
- **Resolving Agents and Solvents chosen on the following criteria:**
 - Literature precedent for optimum crystallisation
 - Compliance/traceability
 - Availability and cost
 - Environmental/FDA regulations (i.e. no acetone due to shipping issues)
- **Ideal for easy robotic manipulations**
 - Kits come with peelable seal or pierceable cap-mat to accommodate direct injection of racemate by liquid dispensers

Types of Screen kits

Primary Screening

- Set of four (4) kits per racemate
 - Comprehensively covers 384 conditions known to offer best separation
 - Racemate/reagent ratio is 1:1
 - Requires totally up to 12 mmol or avg. 3.5 gm racemate (so that diastereomeric salt crystals can be identified with naked eyes; no optical detection needed)
- Kits have 2 components:
 - 96-vial plate has 8 resolving agents (unique reagents in each row); and each vial contains 0.03 mmol of reagent
 - 96-well plate has 12 solvents (unique solvent in each column) which have to be added to the reagent-racemate mixture during reaction

Strong Acid Kits: contain strong acids to be used against weakly base racemate


Reduced Volume Glass Vial Kits: Same as primary screening but totally require only 0.4 mmol or avg. 0.12 gm racemate

- Must be used with HPLC or other optical method of analysis
- Ideal for fully automating screening of library of compounds

How to use ChiroSolv[®] Screen kits

Reaction Flow of Solid Racemate Kit

ChiroSolve, Inc. 



1 Use the ChiroSolv kit for solid racemate




2 Add racemate and remove the "Transfer" solvent



3 Add ChiroSolv solvents



4 Heat the kit until vials contain homogeneous solution



5 Cool the homogeneous content until crystals form



6 Identify vials with crystals for further analysis

Stage-specific chiral resolution needs

Medicinal Chemistry

- Require pure enantiomer for preclinical assessment to avoid confounding effects of testing with racemate
- Fast access to enantiomer with target purity for lead identification

Development Chemistry

- Find efficient, safe, and cost-effective routes for producing enantiomer
- Identify best prospects for scale up and commercial manufacturing
- Adopt 'Green' processes

ChiroSolve medicinal chemistry solution

Medicinal Chemistry

- Require pure enantiomers for preclinical assessment to avoid confounding effects of testing with racemates
- Fast access to enantiomers with target purity for initial testing



- Prepare enantiomer for discovery and preclinical study in matter of days
- Establish early IP foundation by identifying complete *scope of* reagents offering the separation
- Get small quantity of enantiomer for lead identification

Medicinal chemistry: *EnantioPrep*TM

- Prepare pure enantiomer for discovery and preclinical study
- Establish IP scope
- Generate 500 mg to 20 grams of desired enantiomer
- Extract Up to 90% enantiomer from racemate, purity over 95%
- Define the ideal resolution pathway

Separation

Define optimum separation parameters



Recovery

Capture starting material for final separation



Purification

Generate pure enantiomer for testing


Medicinal chemistry: *EnantioPrep*TM

Screen




Matrix of 32 resolving agents and 12 solvents

- Prepare pure enantiomer within days
- Maximize use of racemate
- Establish IP scope
- Generate 500 mg to 20 grams of desired enantiomer
- Define the ideal resolution pathway



Typical: 10 'hits'
Analyze with HPLC



Use the supplied chemicals and tools to recover racemate from all wells



Purify

Use selected resolving agent and solvent and purification tool for enrichment



Pure enantiomer

ChiroSolve development chemistry solution

Development Chemistry


- Find efficient, safe, and cost-effective routes for producing enantiomers
- Identify best prospects for scale up and commercial manufacturing
- Adopt 'Green' processes'



- Efficient routes for producing enantiomer
- Define commercially viable and cost effective processes
- Efficient conversion of existing methods to Green Chemistry processes

Development chemistry - ScalePrep™

- Identify cost-effective and shortest route for manufacturing enantiomers
- Predictive for scaling from kilo lab to ton production
- Optimize process using factors: energy efficiency, safe chemical usage, minimization of waste
- Collect data on alternative processes




Development chemistry: *ScalePrep*TM


Screen



Matrix of 32 resolving agents
and 12 solvents



Typical: 10 'hits';
Analyze with HPLC;
Select best 3 hits



Purify

Using each of 3 combos,
purify and see which offers
pure enantiomer with
minimum steps



Optimize

Order Optimize kits from
Chirosolve to evaluate
optimization factors
(temperature ramps,
stoichiometry)



Optimized crystallization
method

- Choose reagents that offer shortest route to end-point
- Define the ideal resolution conditions for scale-up in terms of proportions of starting material, energy savings

Recovery and Purification kits

After the screening, these kits recover the racemate used; and then obtain pure enantiomer through series of re-crystallization steps. Each kit includes:

- **Recovery solution and recovery unit:** When the contents of the screening kits (after the experiment) is collected into the recovery unit containing the recovery solution, 2 liquid layers are formed: bottom layer contains the racemate; while the top layer contains the resolving agent
- **Purification and Recovery Bottles** that allow easy collection of purified enantiomer and filtrate containing the other enantiomer respectively
- **Filtering unit** with vacuum adapter that allows fast crystal separation during re-crystallization steps

Optimization kits

Based on the results of the screening, customer may request for one or more Optimization Kits. These kits include:

- **Rows of a single resolving agent** in difference quantity to match the ratio 1:1, 1:1.5, 1:2, 1:2.5, 1:3, 1:3.5, 1:4, 1:4.5 with the racemate
- Different columns of the kit containing the **target solvent(s) in different concentration**

Using these kits, process optimization can be achieved within matter of days. These optimization parameters can improve the enantiomeric yield; or can reduce the number of re-crystallization steps. One or more optimization kits containing different combination of reagent/solvents identified as best candidates during screening process can be explored in parallel.


Success stories

- **Discovery group** from **big pharma** was unable to separate out enantiomers; ChiroSolve Screening services were able to give them 4 combinations of resolving agents and solvents that gave good crystals. As a result, they were able to complete the project successfully.
- **Discovery startup** needed a very quick turn-around to explore different separation methods for IP prosecution. Using ChiroSolv kits, they were able to quickly explore the methods without having to compromise quality, leading to patentable results.
- **Process chemistry** group from **major Bio-pharma** used ChiroSolv kits to identify higher purity, and in the process, found a better method that reduced re-crystallization steps from 4 to 3.
- **CRO firm** was able to choose the best derivatisation route that would give optimum end results in terms of enantiomeric purity and yield

Example separations (performed in-house)


- Alcoholic fragment of synthetic pyrethrins**

Allethrin is a potent insecticide and is more stable than the natural analogs. In order to identify the separation process using ChiroSolv, we modified this material with phthalic anhydride & succinic anhydride to prepare acidic monophthalates and succinates.



- n-BOCS-(+)- β -hydroxymethylpiperidine**

This is an intermediate for tryptase inhibitor project.



- Pyridine-containing β -amino acid**

1 (3-pyridyl- β -amino-3-propionic acid) and its ester hydrochloride 2. This combination of acidic and basic functional groups offers a wide choice of resolving agents.

